Validatie

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Verification of qRT-PCR assay for detection of SARS-CoV-2

Method: qPCR – One step RT-PCR with Fast Virus Master Mix

Manufacturer/Origin: Biolegio

Type of validation: verification of improved primers/probe set specific for detection of SARS-CoV-2

Application test: patient diagnostics, detection of SARS-CoV-2 in clinical material

Conclusion: Based on the results, we conclude that using the improved primers/probe set the RT-PCR meets the predefined criteria of the following parameters: measurement trueness, measurement accuracy, quantitation limit, analytical sensitivity and detection limit. The RT-PCR is more sensitive than with the original primers/probe set and now has similar sensitivity as the E-gene RT-PCR. Therefore the improved primers/probe set has been accepted for implementation in routine diagnostics of SARS-CoV-2.


Author: Lisa Wijsman

Date: 6-7-2020

Introduction & purpose

At this moment when a request is made for corona diagnostics the sample is tested on 2 genes of SARS-CoV-2. The E-gene and RdRp-gene are detected in 2 separate PCRs. The E-gene PCR is more sensitive than the RdRp-gene PCR. With an adjustment of the primers and probe set for the RdRp-gene, the idea is that the sensitivity of RdRp-gene PCR improves, so

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Validatie?

Analytic parameters:

The detection qPCR is not used as a pure quantitative method. Results are scored positive/negative to confirm if a patient is infected with SARS-CoV-2. Ct values are used as semi-quantitative information when useful for interpretation of the result. Only those parameters are tested that need verification compared to the original validation

Measurement uncertainty (Meetonzekerheid):

– not tested

Measurement trueness (Juistheid):

– tested: Saliva panel (page 6)

Measurement accuracy (Accuraatheid):

– tested: Saliva panel (page 6)

Measurement precision including measurement repeatability

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Analytical sensitivity (Analytische sensitiviteit):

– tested: LOD95 determination (page 4)

Analytical specificity, including interfering substances

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– tested: 10-fold dilution of quantified control (page 2)

Measuring interval (Meetinterval, kwantitatieve methoden):

– not tested

Clinical parameters:

Diagnostic sensitivity (Diagnostische Sensitiviteit):

– not tested

Diagnostic specificity (Diagnostische Specificiteit):

– not tested